

Protocol for creating unigene set and single seeded decent lines

Plant selection and seed production

Fifty seeds were sown on GM selected media using R7402 20ng/ml and PPT 10µg/ml, after 10-14 days 24 plants per pool were transferred to soil and each plant assigned a unique SM number. Tissue was harvested after 10-14 days and DNA preparations performed using a Qiagen Dneasy 96 plant extraction kit (catalogue number 69181). Seeds were harvested and cleaned after approximately 6 weeks and sent to the Nottingham Arabidopsis Stock Centre (NASC) and assigned a NASC number in addition to a SM number.

Plant selection on soil

Seeds were sown directly on soil and allowed to grow to second leaf stage then sprayed twice at a 7 day interval with BASTA (375mg/ml) and R7402 (1µg/ml or 2µg/ml) and the growth of plants assessed after 7 days.

DNA sequencing and sequence output

DNA samples (2µl) were digested with Apo1 (2U per reaction) for 90 minutes at 50 C and the enzyme inactivated at 75 C for 20 minutes. The samples were ligated overnight (ligase 0.2U at 4 C and 4µl of the ligated DNA aliquoted for the inverse PCR (iPCR) which was performed using two transposon specific primers spm 31 (5'GCT TGT TGA ACC GAC ACT TTT AAC ATA AG 3') and inv6-2 (5' GCT AAG CAC ATA CGT CAG AAA CCA TTA TT 3') at 5 pmoles each and the following cycles, 2 minutes denaturation step at 94 C (94 C 10 seconds 60 C 30 seconds 68 C 5 minutes repeated 38 cycles) and a final step 68 C for 5 minute. iPCR products (20µl) were cleaned up using 2U Exo1 and 2U SAP per reaction at 37 C for 30 minutes and a 10 minute denaturation step at 80 C. The products (5µl) were sequenced using Applied Biosystems Big Dye Terminator kit (catalogue number 4390244) and a nested transposon primer spm 32 (5' TAC GAA TAA GAG CGT CCA TTT TAG AGT GA 3') cycled 25 times at 96 10 seconds, 50 C 5 seconds, 60 C 4minutes. Ethanol precipitated samples were run on a Perkin Elmer 3700 sequencer and results assessed for sequence quality, length (>16 bases), the presence of gene/transposon junction and alignment with Arabidopsis sequence using a p value <0.001.

iPCR protocol

Digestion

DNA	2.0µl (~200ng)
10X Buffer 3 (NEB)	1.0µl
ApoI (4 U/µl)	0.2µl
100XBSA (10mg/ml)	0.1µl
dH ₂ O	6.7µl
TOTAL	10.0µl
Digest at 50°C for 1 hour 30 minutes	
Inactivate enzyme at 75°C for 20 minutes	

Ligation

Digested DNA	10.0µl
10X Ligase Buffer	4.0µl
Ligase (1 U/µl)	0.2µl
dH ₂ O	25.8µl
TOTAL	40.0µl

Ligate overnight at 4°C

PCR

DNA (from ligation)	4.0µl
10X PCR Buffer	2.0µl
dNTPs (2mM)	2.5µl
Primer Spm31 (10µM)	0.5µl
Primer inv6-2 (10µM)	0.5µl
Taq (Qiagen) (1U)	0.2µl
dH ₂ O	10.3µl
TOTAL	20.0µl

Cycling parameters:

1.	94°C	2 min
2.	94 °C	10 sec
3.	60 °C	30 sec
4.	68 °C	5 min
go to 2 for a further 37 cycles		
5.	68 °C	5 min
(extension is 45''-1' per kb)		

Cleanup

PCR products	20µl
SAP (1 U/µl)	2µl
ExoI (10 U/µl)	1µl
Incubate 30 mins at 37°C, 10 mins at 80°C)	

Sequencing reactions

template DNA	5µl
dSpm 32 primer (10µM)	1µl
terminator mix (v3)	3µl
dH ₂ O	1µl

Cycling parameters:

1.	96 °C	10 sec
2.	50 °C	5 sec
3.	60 °C	4 min
go to 2 for a further 24 cycles		

Precipitate sequencing reactions

Add 60µl EtOH (95%) and 1µl EDTA (250mM) to each tube and gently mix
Incubate on ice for 30 mins
Centrifuge at 3500rpm for 30 mins

